



محصولات مراقبت قلبی Teleflex®


تهران، خیابان شهید قندی (پالیزی)، شماره ۱۲۸
کد پستی: ۱۵۵۴۷۴۳۷۶، تلفن: ۸۳۹۶ (خط ۳۰)
No.128, Qandi (Palizi) St., Tehran, Iran
Postal Code:1554744371, Tel: (+9821) 8396
www.iranpharmis.org



Turnpike® LP catheter

RETROGRADE & ANTEGRADE

- دارای انعطاف پذیرترین و سبکترین ساختار در خانواده Turnpike ها
- مناسب برای کیس های Retrograde و Antegrade
- مناسب برای عروق پریفرال و کرونری
- دارای قدرت پیش روندگی بالا به علت داشتن Tip مخروطی و سبک
- دارای شفت ۳ لایه سبک و انعطاف پذیر جهت دسترسی راحت تر به نواحی Distal



[مناسب جهت دسترسی به عروق پر پیچ و غیر مستقیم به علت ساختار بسیار نرم و منعطف]


Navigate Complex Interventions at Every Turn

Ideal Combination of Flexibility and Torque Response

این خانواده از کتترها با طراحی متنوع در ساختار شفت خود رنج وسیعی از انعطاف، چرخش هم محور و استحکام بالا را جهت دسترسی و رفع تصلب های نرم و سخت در عروق کرونری و پریفرال را در اختیار پزشک قرار می دهند.

Outmaneuver Roadblocks

- + Complex lesions
- + Extreme tortuosity
- + Guidewire support



لایه خارجی پلیمری با Taper Tip آبدوست جهت تسهیل حرکت ابزار Intervention

A polymer outer layer, paired with a hydrophilic coating on the distal 60 cm, provides excellent lubricity for delivery

مقاوم در برابر Kink شدن و چرخش غیر هم راستا

Dual-layer coils and a braid make for an ideal combination of flexibility and torque response

حرکت آسان گایدوایر

An inner PTFE liner that spans the entire catheter to ensure smooth guidewire interaction throughout the case


Turnpike® LP

Dual coil tapers down to single coil 21 cm from distal tip for lower profile and increased flexibility in the distal shaft and tip

Ordering Information

MODEL	VERSION	WORKING LENGTH	PROXIMAL SHAFT O.D.	DISTAL SHAFT O.D.	DISTAL TIP O.D.	COMPATIBLE GUIDEWIRE DIAMETER	COMPATIBLE GUIDE CATHETER SIZE
5638	Turnpike LP	135 cm	2.9 Fr. 0.97 mm / 0.038"	2.2 Fr. 0.74 mm / 0.029"	1.6 Fr. 0.53 mm / 0.021"	0.014" (0.36 mm)	5 Fr. 1.42 mm / 0.056" I.D.
5639	Turnpike LP	150 cm	2.9 Fr. 0.97 mm / 0.038"	2.2 Fr. 0.74 mm / 0.029"	1.6 Fr. 0.53 mm / 0.021"	0.014" (0.36 mm)	5 Fr. 1.42 mm / 0.056" I.D.

Packaged in quantities of 1 unit per box.



1 All values based on bench test data averages compared to the standard Turnpike Catheter, n=15, performed by Teleflex. Bench test results may not necessarily be indicative of clinical performance. Data on file.

2 All values based on bench test data averages, n=2, performed by Teleflex. Bench test results may not necessarily be indicative of clinical performance. Data on file.